

With Swinging instruments hoses :

It's a new generation unit, a combination of the high Swiss technology with Italian production and design. Ergonomic handling, smooth operating and a high quality at very competitive price.



The Electromechanical Patient Chair "AG03SWF" with motor lifting to low tension

- The headrest is provided with a double articulation
- 4=3 different chair programs positions and plus one additional rising and return to the last work position
- All positions are stored and controlled from the dentist side and part from the assistant's side
- The chair is equipped with 2 safety switches
- No cables or supply hoses on floor

The Hydro Unit :

- Removable and auto-cleavable ceramic cuspidor bowl
- Automatic timer of cup filler and bowl flush.
- Tube of 20" for a direct discharge
- Filters for aspiration and cuspidor bowl are inspected from the outside part

Spécifications Techniques/Détails :

Dental Unit Friend Up Art .Nr. SW77120 with Swinging hoses

system: complete with (Instrumentation From left to right):

- Pos.1=1 3-way f syringe,
- Pos.2=Turbine tubing with light
- Pos.3=Micro motor tubing with micro motor MC2 LK ,
- Pos.4=Free (2nd Turbine tubing or Micro motor optional)
- Pos.5=Free (EMS or SUPPRASSON Scaler or curing light Mini Led optional)
- Control panel for patient's chair and instruments
- Swivel tray support with one tray

Hydro Group unit :

- Removable and Auto-cleavable ceramic cuspidor bowl
- Can rotate 180° (space saving)
- Cuspidor's drainage valve included
- Tube of 20" for a direct discharge
- **Safe Spray** distilled water system (Optional)

Assistant table : (Instruments from left to right)

- Extensions options to 3 instruments
- Pos.1=1 Surgical section,
- Pos.2=1 saliva ejector,
- Pos.3=Free (optional= 3-way-syringe for assistant or curing light)
- **FS40** The control panel with 2 Program Buttons for Flush bowl and cup filler

Fauteuil patient AG03SWF:

- 4=3 programs positions, plus additional rising position and return to the last position.
- The Head rest with a double articulation

The foot control:

- The foot control for the command of Instruments functions.
- Activation of the instruments with a progressive regulation of micro motors power.
- The spray and chip blower function are included.

The operating Lamp: (Unit Mounted)

- Equipped with EDI Light (MAIA or ALYA LED light optional)
- Reliable and proved
- 24.000 Lux, 5000 °k

Colors :

- For exact definition colors please refer to the separate Chair's upholstery color table
- Unit base, arm systems, and Hydro unit standard color :(color choose optional)

